Methods for leaf-level photochemistry and gas exchange

Copyright © 2024 LI-COR, Inc. All rights reserved. Contact LI-COR for permission to redistribute.



Some opening remarks

- Much of this talk is organized around assumption and potential errors in measurements
- Review of photochemistry and why one might want to measure it
- Deep dive into instrument theory behind the techniques
- Some practical considerations for making measurements
- Demonstration measurement



































How do we measure photochemistry?





Combined gas exchange and fluorescence instruments

LI-600/N Porometer Fluorometer

> LI-6800 Portable Photosynthesis System

The LI-6400/XT is still widely used, but we don't manufacture it anymore. So, it doesn't get to be in the picture!



Some physiological parameters from gas exchange

Outputs from a gas exchange system









Parameters mean different things depending on how the instrument is applied

Survey: Tells us about the current state of the pant







Example response to light (AQ curve)





Example response to light (AQ curve)





Example response to light (AQ curve)





Example response to CO_2 (AC_i curve)





Example response to CO₂ (AC_i curve)





Why measure photochemistry and gas exchange?

The processes these measurements probe are processes fundamental to support (almost) all life on Earth



Traditional breeding and selection





Traditional breeding and selection





The plant as a factory





The plant as a factory





How do we measure leaf-level photochemistry and gas exchange?



Energy capture and dissipation

Idealized Jablonski diagram of states and dissipation at the single molecule level





Fluorescence yield as defined by dissipation path



 f_F = Flux of fluoresced photons

 k_x = Rate constant of dissipation path

q = Absorbed photons Φ_F = Fluorescence yield

$$f_F = q \frac{k_F}{k_D + k_F + k'_{PC} + k'_{NPQ}}$$

$$\phi_F = \frac{f_F}{q}$$

$$\phi_F = \frac{k_F}{k_D + k_F + k'_{PC} + k'_{NPQ}}$$

28 LI-COR.

Fluorescence yield as defined by dissipation path



$$\phi_F = \frac{k_F}{k_D + k_F + k'_{PC}p_{Q_A} + k'_{NPQ}p_z}$$

Photochemistry and NPQ are apparent rate constants, adjusted by proportionality coefficients:

 p_z is the proportion of NPQ engagement

 p_{Q_A} is related to the proportion of open reaction centers

LI-COR.

29

 f_F = Flux of fluoresced photons q = Absorbed photons Φ_F = Fluorescence yield k_x = Rate constant of dissipation path

Pulse Amplitude Modulation (PAM) fluorometry











Time (not to scale)



Deriving fluorescence parameters

$$\phi_F = \frac{k_F}{k_D + k_F + k'_{PC}p_{Q_A} + k'_{NPQ}p_z}$$

There are assumptions of state for each parameter that can be defined using the proportionality coefficients. Substituting these definitions for each parameter...



Deriving fluorescence parameters

$$F_o = \frac{k_F}{k_D + k_F + k'_{PC}}$$

$$F_m = \frac{k_F}{k_D + k_F}$$

$$F_s = \frac{k_F}{k_D + k_F + k'_{PC}p_{QA} + k'_{NPQ}p_z}$$

$$F'_m = \frac{k_F}{k_D + k_F + k'_{NPQ}p_z}$$





Key assumptions





Modulation for F_0 (verify $p_{Q_A}=1$)

The higher the modulation frequency the larger the signal and the smaller the noise.




Test for saturating flash intensity (verify $p_{Q_A}=0$)



See Markgaf T and J Berry. 1990. Measurement of photochemical and non-photochemical quenching: Correction for turnover of PS2 during steady-state photosynthesis. Research in Photosynthesis

Multi-Phase Flash (MPF) technique

- Segment a sub-second flash into a three phase protocol
- Phase 1 and 3 look like ends of a traditional square flash
- Phase 2 includes a linear ramp in intensity
- Fit phase 2 fluorescence versus the inverse of intensity to estimate F'_m at infinite light intensity

Loriaux SD *et al.* 2013. Closing in on maximum yield of chlorophyll fluorescence using a single multiphase flash of sub-saturating intensity. Pant, Cell and Environment



40 LI-COR.

Methodologies for determining F'_o

Methodologies for determining F'_o

$$F_o' = \frac{k_F}{k_F + k_H + k_{PC}' + k_{NPQ}' p_Z}$$

NPQ has both slow and fast relaxing components. In some plants the fast components are fast enough to invalidate this assumption. Where this happens, it will be particularly evident at high light where NPQ makes up a large portion of the quenching capacity.

NPQ is left intact: 0<pz≤1

Oxborough and Baker. 1997. Resolving chlorophyll a fluorescence images of photosynthetic efficiency into photochemical and nonphotochemical components – calculation of qP and Fv-/Fm-; without measuring Fo-. Photosynthesis Research

Fluorescence is an uncalibrated measurement

- Meaningful parameters are all ratios
- Signal is a raw count, and not standardized to any specific scalar
- Only a fraction of the total emission is seen
- Emission characteristics are dependent on sample optical properties

Light absorption and weighted absorption coefficients

- Absorption is measured with a spectroradiometer and integrating sphere (black curve)
- Absorption coefficients weighted by the spectral output of the light source in particular bands (α_x)
- Absorbed light accounts for the fraction of light output in particular bands (p_x)

$$Q_{absorbed} = Q_{incident}(\alpha_{blue}p_{blue} + \alpha_{red}p_{red})$$

 Coefficients are calculated by integrating the product of the normalized output for each band with the absorption curve

Data from an algal suspension at a typical 6800-18 measurement density. LED normalized output for 6800-01A.

Contact envsupport@licor.com and provide light source model number for spectral outputs.

8/18/2024

44

LIEDR

What's what in a LI-COR gas exchange system

- Sensor head contains the core of the gas exchange sensors
- Console manages the air supply, contains the power supply and includes the user interface

In the LI-600 the "sensor head" and "console" are effectively merged into one package

What's what in a LI-COR gas exchange system

- Sensor head contains the core of the gas exchange sensors
- Console manages the air supply, contains the power supply and includes the user interface

In the LI-600 the "sensor head" and "console" are effectively merged into one package

Gas exchange parameters

Outputs from a gas exchange system

s = Leaf surface area $f_x = \text{Flux of gas } \times$ $u_{in} = \text{Molar air flow entering the chamber}$ $u_{out} = \text{Molar air flow exiting the chamber}$ $x_{in} = \text{Mole fraction of } \times \text{ entering chamber}$ $x_{out} = \text{Mole fraction of } \times \text{ exiting the chamber}$

$$sf_x = u_{out}x_{out} - u_{in}x_{in}$$

This is the model that has been used in the LI-6400/XT and LI-6800 traditionally and what the LI-600 uses

Steady state versus dynamic

Example data from an empty chamber (Assimilation = 0)

At steady-state - Adyn 60 - Asty 520<u>CO2_r</u> 520<u>CO2_s</u> the transient term is zero and 40 500 the fluxes match 20 480 S^{-1} µmol molµm9/m-2 460 -20 440 -40 420 -60 400 -2 -1.5 -0.5 -1 Meas:TIME(min)

Following a step change in CO_2 the steady-state model gives a false flux driven by the volumetric time constant.

The dynamic model accounts for this and still shows zero assimilation.

Steady state versus dynamic

- Biology versus instrument
- Characterize fast biological transient responses
- Non-steady state response curves
 - High data density in less time

Fluxes are diffusive in nature

Carbon dioxide and water vapor move by diffusion between the well mixed atmosphere and the inside of the leaf, each driven by a concentration gradient:

 $f_x = -D_x \frac{\partial x}{\partial d} = g_{t_x}(x_1 - x_0)$

Fluxes are diffusive in nature

$$E \approx g_{t_w}(w_i - w_a)$$
$$A \approx g_{t_c}(c_a - c_i)$$

The approximation here is because these ignore diffusive collisions with water and air, and the transpiratory flux impact on assimilation. 8/18/2024

Fluxes are diffusive in nature

The approximation here is because these ignore diffusive collisions with water and air, and the transpiratory flux impact on assimilation.

e

$$w_{i} = \frac{e(T_{leaf})}{P} 1000$$
Neasured
fluxes

$$W_{i} = \frac{W_{out}}{P} 1000$$

$$W_{a} = W_{out}$$

$$W_{a} = W_{out}$$

$$(E) \approx g_{t_{w}} (W_{i}) - W_{a}$$

$$(G_{a}) = C_{i}$$

$$g_{t_{c}} = \frac{g_{t_{w}}}{1.6}$$

$$c_{a} = C_{out}$$

Separating stomatal contribution

The diffusive pathway can be divided into segments, each making some contribution to the total conductance. Summing those contributions follows Ohm's Law:

$$g_t = \frac{1}{r_t} = \frac{1}{r_{bl} + \left(\frac{1}{r_s + r_i} + \frac{1}{r_c}\right)^{-1}}$$

To get the stomatal contribution, some simplifying assumptions are made:

- Intercellular resistance is small
- Cuticular resistance is large

$$g_{s_w} = \frac{1}{\frac{1}{g_{t_w}} - \frac{1}{g_{bl_w}}}$$

Separating stomatal contribution

The discussion so far, only considers a single leaf surface, analogous to how conductance is measured in the LI-600. For two leaf surfaces, as measured in the LI-6800 or LI-6400/XT, we account for the effect of stomata on both surfaces with an adjustment to the boundary layer contribution:

Boundary layer conditions

Leaf in the bulk atmosphere

Boundary layer conductance is inversely proportional to its thickness.

Despite typically high(er) ambient VPD (or low relative humidity), leaves in the bulk atmosphere generally experience something relatively close to saturation at the leaf surface.

Boundary layer conditions

Leaf in the bulk atmosphere

Leaf in the chamber

What does this mean for making measurements?

Target boundary layer conditions!

$$RH = \frac{e}{e_{T_{air}}} 100$$

$$VPD_{air} = e_{T_{air}} - e$$

$$VPD_{leaf} = e_{T_{leaf}} - e$$

On the LI-6400/XT when using automated control, W_{out} (H2OS) will give the best control. On the LI-6800 use constant VPD_{leaf} (1 to 1.5 kPa)

RH = Relative humidity

e = Vapor pressure

 e_{T_x} = Saturation vapor pressure at temperature T_x

 \tilde{VPD}_x = Vapor pressure deficit

The apparent nature of transpiration?

The top panel is modeled assuming a constant conductance and ignoring boundary layer effects. The latter only shifts the curves left or right, it doesn't change their shape.

The bottom panel is actual leaf data collected by varying the chamber humidity across a range of "typical" boundary layer humidity.

What should we conclude from this?

Transpiration is a real count of water vapor molecules leaving the leaf.

Transpiration is artifactual in nature, driven by the chamber environment as much as by biology.

Conductance removes the environmental dependence

_

What does this mean for interpreting the data?

Be careful trying to draw conclusions based on apparent transpiration! Water use efficiencies (WUE) provide a nice illustration of this...

$$W_t = \frac{A}{E} = \frac{A}{g_t(w_i - w_a)}$$

$$g_{s_w}$$

 $W_{a} =$

For a nice overview of WUE parameters see Seibt et al. 2008. Carbon isotopes and water use efficiency: sense and sensitivity. Oecologia

 W_t = Instantaneous WUE W_a = Intrinsic WUE

What about *actual* (evapo)transpiration?

There are other, less constrained ways to get to this number...

- Eddy covariance or similar methods: Allow for continuous monitoring
- Permits undisturbed boundary layer conditions
- Field/ecosystem scale measurements

Some practical considerations for measurements

The leaf's light environment

Let's consider this plant canopy and light reaching it... Note the clearly demarcated shadows.

Page 67

The leaf's light environment

LI-190/R

- Filtered to idealized plant spectral response
- Measures light incident on the plane of the diffusor

Page 68

Controlling CO₂ at the leaf

You can either keep the incoming or the out-going (chamber) CO₂ concentration constant, but not both. Whichever one is not held constant will vary proportionally with the assimilation rate. So, sensitivity to CO₂ concentration is an important consideration!

Flow control: optimize for the "signal"

Fixed flow rate

Choose the highest flow rate that still gives a good delta

Humidity control

Hold VPD_{leaf} between 1 to 1.5 kPa

Data from an LI-6400XT and 6400-02B chamber. Error bars are an estimated uncertainly assuming typical noise on the LI-6400 IRGAs.

Time of day effects

At some point in the day, the correlation between assimilation and PAR breaks. This typically happens near or after peak light intensity is reached and is referred to as the mid-day depression.

Selecting an appropriate chamber?

What specific data do you need?

What about leaves that don't fill the chamber?

Where the leaf doesn't fill the chamber, area needs to be measured and the correct leaf area used in the calculations

In the fluorometer, the bottom of the chamber is removable allowing the leaf to be imaged while its still in the chamber





73

What about leaves that don't fill the chamber?

The LI-600 chamber must be filled!







74

Not just leaves!



6800-18 Aquatic Chamber

For measurements on liquid samples or water saturated materials



Questions?

This document is really useful! Includes a description of parameters with relevant literature references for each. https://licor.app.box.com/s/wcglihmvd1rwotm0i0avh5vwd0r4lb7g

Complementary Leaf Physiology Measurements: The LI-600 and LI-6800

Ø Application Note

Contents

Abstract
Gas exchange and photosynthesis
The LI-800 Parameter/Fluorometer
The LI-6800 Photosynthesis System
System comparison
Parameters measured
Case study 1: Water use in wild peppers
Case study 2: Parameterizing stomatal ratio

Abstract

The LI-800 Porometer/Fluorometer is designed for highspeed sampling. It quickly measures stomatal conductance and chlorophyll a fluorescence from leaves. In contrast, the LI-8800 Portable Photosynthesis System is designed to provide highly detailed data. It measures parameters that take more time than porometry alone, such as $\rm CO_2$ assimilation. Although the LI-800 and LI-8800 have distinct purposes, the data they provide is highly complementary. When used together, they can improve the efficiency of data collection in laboratories, greenhouses, and plant breeding programs.

In this document, we discuss how these instruments can be used together to collect larger data sets and explore physiological traits in more detail while taking less time. To illustrate the concept, we review a case study in which wild pepper varieties were screened for sensitivity to water stress using the LI-800 and measured in detail with the LI-8800. Over 500 individual leaves were measured in a single day using one LI-800. Data from the high-volume survey were used to identify individuals that were evaluated in detail using the LI-8800.

Gas exchange and photosynthesis

Leaves perform many complex physiological processes, but often of particular interest is how they exchange carbon dioxide (CO₂) and water (H₂O) between the leaf interior and the atmosphere (summarized in Figure 1). These represent the two major processes, assimilation (CO₂) and transpiration (H₂O), and the movement of these gases into and out of the leaf is typically facilitated by the stomates. The stomates are able to open or close in response to stimuli or stress in order to control the rate of exchange.

Along with the diffusion of gases into and out of the leaf, we have the photosynthetic apparatus within the chloroplasts. Here, energy from captured photons facilitates the conjugation of CO_2 into sugars. This process consists of the light-independent reactions or Calvin-Benson cycle occurring in the stroma and the light-dependant reactions occurring in the thylakoid membranes.





Figure 1. Cross section of a leaf (top) and cell (bottom) showing a summary of the biochemical processes within. These processes are probed by the LI-600 and LI-6800.

76

LI-COR.

